



PRODUCT COMPARISON

FEATURES	ARTEL PIPETTE TRACKER™	ARTEL PIPETTE TRACKER PRO™
Bar code scanning to quickly and accurately identify pipettes*	Yes	Yes
Calibration labels and bar code identity labels*	Yes	Yes
Trending reports for pipette calibration performance over time	Yes	Yes
Default databases in Microsoft® Access	Yes	Yes
Manual data entry	Yes	Yes
Redo sample option	Yes	Yes
Records date of pipette entry	Yes	Yes
Imports pipette entries from text files	Yes	Yes
Includes Artel tolerance limit table	Yes	Yes
Preview of system and calibration reports	Yes	Yes
Filterable data selection	Yes	Yes
Export data in HTML, PDF, XLS, CSV, TIF, RTF formats	Yes	Yes
Email notification when pipette calibration is due	No	Yes
Audit trail system & ability to network databases	No	Yes
Tracks and saves updated protocols	No	Yes
Facilitates 21 CFR Part 11 compliance for security features	No	Yes
Administrative review and approval of calibration reports with ability to enter comments post calibration.	No	Yes
Supports SQL or ORACLE® databases	No	Yes

* Requires the purchase of PCS-795: Pipette Tracker Enhancements Hardware Package including Label Printer and Bar Code Reader.