



ARTEL Pipette Tracker™

Comparison

FEATURES	ARTEL PIPETTE TRACKER™	ARTEL PIPETTE TRACKER™ PRO
Bar code scanning to identify pipettes and launch a calibration*	✓	✓
Print bar code identity and calibration labels displaying Calibrated by, Status and Next due date*	✓	✓
Prevent As-Found failures with Action Limits to alert you if results are approaching tolerance limits	✓	✓
Maintenance plan to record and report service & repairs	✓	✓
Trending reports for pipette calibration performance over time	✓	✓
Databases in Microsoft® SQL Server Express 2005 or ORACLE®	✓	✓
Manual data entry for all volumes	✓	✓
Option to redo or delete sample	✓	✓
Records date of pipette entry	✓	✓
Imports pipette entries from text files	✓	✓
Includes Artel recommended tolerance limits	✓	✓
Preview of system and calibration reports	✓	✓
Filterable data selection	✓	✓
Export data in HTML, PDF, XLS, CSV, TIF, RTF formats	✓	✓
Email notification when pipette calibration is due		✓
Audit trail system & ability to network databases		✓
Tracks and saves updated protocols		✓
Facilitates 21 CFR Part 11 compliance for security features		✓
Administrative review and approval of calibration reports with ability to enter comments post-calibration		✓

* Requires the purchase of PCS-795: Pipette Tracker Enhancements Hardware Package including Label Printer and Bar Code Reader.

19A5117F

25 Bradley Drive, Westbrook ME 04092-2013
 toll-free: 888-406-3463 • tel: 207-854-0860 • fax: 207-854-0867 • email: info@artel-usa.com

www.artel-usa.com