

ARTEL Pipette Calibration Report

Pipette ID: 100-1002
Description: Gilson Pipetman 100, variable, 100µL, 1 channel

PCS Instrument Details

User Entry

Reagent Lot Code: 12345
Temperature: 22.5 °C

Location: ARTEL Laboratory
Owner:
Method: Control Calibration
Test Plan: 100 uL, fixed
Type: As Calibrated
Interval: 02
Comment:

Pipette Tip Information

Tip Name: Gilson 100µL
Tip Lot: 3106 ALTS OV

STATUS: Passed

Target Volume (µL)	Channel	Mean Volume (µL)	Standard Deviation	Relative Inaccuracy (%)		CV(%)		Status
				Actual	Target	Actual	Target	
100.00	1	99.8000	0.1370	-0.20	2.00	0.14	1.00	Passed

Sample	Volume (µL)
1	99.9800
2	99.8900
3	99.6300
4	99.6500
5	99.7600
6	99.8900
7	99.9500
8	99.7800
9	99.8400
10	99.5900

Performed by: ARTEL Technician

Date: 08 February 2008 at 10:57:26

Reviewed by: _____

Date: _____

Printed Date: 2/8/2008

Printed Time: 10:57:34 AM